


# Homocysteine

Enzyme Cycling Sequence  
for Greater Sensitivity and  
Stability



Calibrator Standardized to  
NIST SRM 1955 "Homocysteine  
Standards Reference Material"

- Excellent Correlation to HPLC & EIA
- EDTA or SST Samples Acceptable
- Extended / Optimized Calibration Stability
- Up to 30 day Reagent On-board Stability
- Two (2) Reagents: Liquid, Ready-to-Use
- Sensitivity 1 $\mu$ mol/L
- Two (2) Calibrators: Liquid, Ready-to-Use
- Controls: Liquid, Ready-to-Use
- Instrument Applications Available

## ORDERING INFORMATION

H7575-57: Homocysteine  
Reagent Set  
With Calibrators  
R1: 1x50ml  
R2: 1x7ml  
Calibrators: 2x3ml

H7575-CTL: Homocysteine Control Set  
(Tri-Level)  
L1: 2X5ml  
L2: 2X5ml  
L3: 2x5ml

CPT CODE: 83090  
Reimbursement: ~ \$23.50 / U.S. Only

*Reliable Reagents • Reliable Results*

For more Information contact:  
Larry Noel, U.S. Sales Manager  
Sales@pointscientific.com  
www.pointscientific.com  
800-445-9853

For International Sales please contact:  
Cathie Kowalski, International Sales  
ckowalski@pointscientific.com  
www.pointscientific.com