

PRESS RELEASE

DATE: September 16, 2009
CONTACT: Ms. Cathie Kowalski, International Sales
PHONE: 734-487-8313
EMAIL: ckowalski@pointescientific.com
WEB: www.pointescientific.com

FOR IMMEDIATE RELEASE

New CK-MB Formulation now available

Pointe Scientific, Inc. has introduced a new formulation of CK-MB reagent. CK-MB activity increases significantly following acute myocardial infarction. This elevation is highly specific for the laboratory diagnosis of myocardial infarction, commonly known as heart attack. Blood tests for CK-MB are commonly used to help confirm myocardial infarction or other distresses to the heart.

Pointe Scientific's new CK-MB formulation does not require the use of a diluent, and is reconstituted with distilled water. The elimination of the diluent will decrease shipping weight and cost, this is idea where transport costs are relatively higher. The reagent is easier to use, and is stable for 5 days at 2-8°C once reconstituted. Catalog numbers include C7562-26, C7562-65, and C7562-150. Controls are also available. For more information please contact your sales representative or our Technical department via email at tech@pointescientific.com or by phone at 734-487-8319.

#####