



Albumin

METHOD

Bromcresol Green, Endpoint

TECHNICAL INFORMATION

Storage Temperature:
Reconstituted Stability
Linearity:
Wavelength:
Expected Values

RT
Ready-to-use Liquid
0.5-8.0 g/dL
630 nm
3.5-5.3 g/dL

ORDERING INFORMATION

CATALOG #
A7502-120
A7502-500
A7502-1L
Calibrator/Standard Required

SIZE
1 x 120 ml
1 x 500 ml
1 x 1 L

RELATED PRODUCTS

A7502-STD
C7506-50

Albumin Standard
Chemistry Calibrator