



Alcohol

METHOD

Alcohol Dehydrogenase, Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	7 Days 2-8°C
Linearity:	0.4 g/dL
Wavelength:	340 nm
Expected Values:	<0.1g/dL

ORDERING INFORMATION

CATALOG #	SIZE
A7504-39	6 x 6.5 ml
A7504-150	10 x 15 ml
Standard Required	

RELATED PRODUCTS

A7504-STD	Alcohol Standard
A7504-CTL	Ammonia/Alcohol Control Set