



Alkaline Phosphatase

METHOD

p-NPP/Kinetic

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	7 Days RT 60 Days 2-8°C
Linearity:	800 U/L
Wavelength:	405 nm
Expected Values:	Adults 35-123 U/L
Children have a higher normal value	

ORDERING INFORMATION

CATALOG #	SIZE
A7505-150	10 x 15 ml
A7505-300	20 x 15 ml
A7505-400	8 x 50 ml
A7505-1L	4 x 250 ml