



## Alkaline Phosphatase Liquid

### METHOD

p-NPP, Kinetic

### TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability	Ready-to-use Liquid
Linearity:	0-1000 U/L
Wavelength:	405 nm
Expected Values:	Adults 35-123 U/L
Children have a higher normal value	

### ORDERING INFORMATION

CATALOG #	SIZE
A7516-150	1 x 120 ml (R1) 1 x 30 ml (R2)
A7516-450	3 x 120 ml (R1) 1 x 90 ml (R2)
A7516-625	1 x 500 ml (R1) 1 x 125 ml (R2)