



ALT (SGPT)

METHOD

Optimized Tris Buffer, Kinetic

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	48 Hours RT 30 Days 2-8°C
Linearity:	500 U/L
Wavelength:	340 nm
Expected Values:	Up to 38 U/L

ORDERING INFORMATION

CATALOG #	SIZE
A7525-75	5 x 15 ml
A7525-150	10 x 15 ml
A7525-300	20 x 15 ml
A7525-400	8 x 50 ml
A7525-1L	4 x 250 ml