



ALT (SGPT) Liquid

METHOD

Modified IFCC, Kinetic

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	0-500 U/L
Wavelength:	340 nm
Expected Values:	4 to 36 U/L

ORDERING INFORMATION

CATALOG #	SIZE
A7526-150	1 x 120 ml (R1) 1 x 30 ml (R2)
A7526-450	3 x 120 ml (R1) 1 x 90 ml (R2)
A7526-625	1 x 500 ml (R1) 1 x 125 ml (R2)