



Ammonia

METHOD

∞ Ketoglutarate, Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	15 Days 2-8°C
Linearity:	600 umol/L
Wavelength:	340 nm
Expected Values:	11-35 umol/L

ORDERING INFORMATION

CATALOG #		SIZE
A7553-34	Substrate	4 x 6.5 ml
	Enzyme	4 x 2 ml
A7553-85	Substrate	10 x 6.5 ml
	Enzyme	10 x 2 ml

Standard Required

RELATED PRODUCTS

A7553-STD100	Ammonia Standard
A7553-STD500	Ammonia Standard
A7504-CTL	Ammonia/Alcohol Control Set