



**AST (SGOT) Liquid**

**METHOD**

Modified IFCC, Kinetic

**TECHNICAL INFORMATION**

|                          |                     |
|--------------------------|---------------------|
| Storage Temperature:     | 2-8°C               |
| Reconstituted Stability: | Ready-to-use Liquid |
| Linearity:               | 0-500 U/L           |
| Wavelength:              | 340 nm              |
| Expected Values:         | 5-34 U/L            |

**ORDERING INFORMATION**

| CATALOG # | SIZE                               |
|-----------|------------------------------------|
| A7561-150 | 1 x 120 ml (R1)<br>1 x 30 ml (R2)  |
| A7561-450 | 3 x 120 ml (R1)<br>1 x 90 ml (R2)  |
| A7561-625 | 1 x 500 ml (R1)<br>1 x 125 ml (R2) |