



Amylase

METHOD

CNPG3, Kinetic

TECHNICAL INFORMATION

Storage Temperature:
 Reconstituted Stability:
 Linearity:
 Wavelength:
 Expected Values:

2-8°C
 Ready-to-use Liquid
 0-2,000 U/L
 405 nm
 25-125 U/L

ORDERING INFORMATION

CATALOG #
 A7564-60
 A7564-120
 A7564-480
 A7564-1L

SIZE
 2 x 30 ml
 4 x 30 ml
 4 x 120 ml
 1 x 1 L