



## Bilirubin, Direct

### METHOD

Diazo,Endpoint

### TECHNICAL INFORMATION

|                          |                     |
|--------------------------|---------------------|
| Storage Temperature:     | RT                  |
| Reconstituted Stability: | Ready-to-use Liquid |
| Linearity:               | 20 mg/dL            |
| Wavelength:              | 550 nm              |
| Expected Values:         | 0-0.5 mg/dL         |

### ORDERING INFORMATION

| CATALOG # |         | SIZE       |
|-----------|---------|------------|
| B7538-120 | Reagent | 1 x 120 ml |
|           | Nitrite | 1 x 15 ml  |
| B7538-500 | Reagent | 1 x 500 ml |
|           | Nitrite | 1 x 15 ml  |
| B7538-1L  | Reagent | 1 x 1 L    |
|           | Nitrite | 1 x 15 ml  |

Calibrator Required

### RELATED PRODUCTS

|           |                      |
|-----------|----------------------|
| B7538-CAL | Bilirubin Calibrator |
| C7506-50  | Chemistry Calibrator |