



BUN (Berthelot) Colorimetric

METHOD

Berthelot, Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	70 mg/dL
Wavelength:	630 nm
Expected Values:	7-18 mg/dL

ORDERING INFORMATION

CATALOG #	SIZE
B7551-120	120 test
Enzyme	1 x 60 ml
Color	1 x 60 ml
Base	1 x 240 ml
Standard	1 x 4 ml

RELATED PRODUCTS

B7550-STD BUN Standard