



BUN (Urea Nitrogen) Liquid

METHOD

Urease/GLD, Initial Rate

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	0-150 mg/dL
Wavelength:	340 nm
Expected Values:	7-18 mg/dL

ORDERING INFORMATION

CATALOG #	SIZE
B7552-150	1 x 120 ml (R1) 1 x 30 ml (R2)
B7552-450	3 x 120 ml (R1) 1 x 90 ml (R2)
B7552-625	1 x 500 ml (R1) 1 x 125 ml (R2)

Calibrator/Standard Required

RELATED PRODUCTS

B7550-STD	BUN Standard
C7506-50	Chemistry Calibrator