



Bilirubin, Total

METHOD

Walters & Gerarde, Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	20 mg/dL
Wavelength:	550 nm
Expected Values:	0.2-1.2 mg/dL

ORDERING INFORMATION

CATALOG #		SIZE
B7576-120	Reagent	1 x 120 ml
	Nitrite	1 x 15 ml
B7576-500	Reagent	1 x 500 ml
	Nitrite	1 x 15 ml
B7576-1L	Reagent	1 x 1 L
	Nitrite	1 x 15 ml

Calibrator Required

RELATED PRODUCTS

B7538-CAL	Bilirubin Calibrator
C7506-50	Chemistry Calibrator