



# Chloride

## METHOD

Mercuric Thiocyanate, Endpoint

## TECHNICAL INFORMATION

Storage Temperature:	RT
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	80-120 mEq/L
Wavelength:	480 nm
Expected Values:	98-106 mEq/L

## ORDERING INFORMATION

CATALOG #		SIZE
C7501-120	Reagent	1 x 120 ml
	Blank	1 x 15 ml
C7501-500	Reagent	1 x 500 ml
	Blank	1 x 30 ml
C7501-1L	Reagent	1 x 1 L
	Blank	2 x 30 ml

*Calibrator/Standard Required*

## RELATED PRODUCTS

C7501-STD	Chloride Standard
C7506-50	Chemistry Calibrator