



Calcium (2 part liquid)

METHOD

Cresolphthalein Complexone, Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	20 mg/dL
Wavelength:	570 nm
Expected Values:	8.5-10.5 mg/dL

ORDERING INFORMATION

CATALOG #		SIZE
C7503-120	Color	2 x 30 ml
	Buffer	2 x 30 ml
C7503-480	Color	2 x 120 ml
	Buffer	2 x 120 ml

Calibrator/Standard Required

RELATED PRODUCTS

C7503-STD	Calcium Standard
C7506-50	Chemistry Calibrator