



## Calcium (Dry)

### METHOD

Cresolphthalein Complexone, Endpoint

### TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	5 Days RT 30 Days 2-8°C
Linearity:	20 mg/dL
Wavelength:	570 nm
Expected Values:	8.5-10.4 mg/dL

### ORDERING INFORMATION

CATALOG #	SIZE
C7508-300	20 x 15 ml
C7508-400	8 x 50 ml
<i>Calibrator/Standard Required</i>	

### RELATED PRODUCTS

C7503-STD	Calcium Standard
C7506-50	Chemistry Calibrator