



Creatine Kinase (CK)

METHOD

CK-NAC, Kinetic

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	24 Hours RT 30 Days 2-8°C
Linearity:	1,500 U/L
Wavelength:	340 nm
Expected Values:	
Males:	Up to 160 U/L
Females:	Up to 130 U/L

ORDERING INFORMATION

CATALOG #	SIZE
C7512-39	6 x 6.5 ml
C7512-120	20 x 6.5 ml
C7512-300	20 x 15 ml
C7512-400	8 x 50 ml
C7512-1L	4 x 250 ml