



Creatine Kinase (CK) Liquid C7522

Method CK-NAC, Kinetic

Technical Information

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	1,200 U/L
Wavelength:	340 nm
Expected Values:	
Males:	38 - 174 U/L (37°C)
Females:	26 - 140 U/L (37°C)

Ordering Information

CATALOG#	SIZE
C7522-150	1 x 120 ml (R1) 1 x 30 ml (R2)
C7522-450	3 x 120 ml (R1) 1 x 90 ml (R2)
C7522-625	1 x 500 ml (R1) 1 x 125 ml (R2)