



## Calcium (Arsenazo III)

### METHOD

Arsenazo III, Endpoint

### TECHNICAL INFORMATION

Storage Temperature:	RT
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	0-15 mg/dL
Wavelength:	650 nm
Expected Values:	8.5-10.4 mg/dL

### ORDERING INFORMATION

CATALOG #	SIZE
C7529-120	1 x 120 ml
C7529-500	1 x 500 ml
C7529-1L	1 x 1 L
<i>Calibrator/Standard Required</i>	

### RELATED PRODUCTS

C7503-STD	Calcium Standard
C7506-50	Chemistry Calibrator