



Creatinine

METHOD

Modified Jaffe, Initial Rate

TECHNICAL INFORMATION

Storage Temperature:	RT
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	25 mg/dL
Wavelength:	510 nm
Expected Values:	0.40-1.40 mg/dL

ORDERING INFORMATION

CATALOG #	SIZE
C7539-150	1 x 120 ml (R1) 1 x 30 ml (R2)
C7539-625	1 x 500 ml (R1) 1 x 125 ml (R2)
C7539-1250	1 x 1 L ml (R1) 1 x 250 ml (R2)

RELATED PRODUCTS

C7513-STD	Creatinine Standard
C7506-50	Chemistry Calibrator