



Creatine Kinase - MB

METHOD

Immunoinhibition, Kinetic

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	5 Days 2-8°C
Linearity:	≤175 U/L
Wavelength:	340 nm
Expected Values:	0-24 U/L

ORDERING INFORMATION

CATALOG #		SIZE
C7562-26	Reagent	4 x 6.5 ml
C7562-65	Reagent	10 x 6.5 ml
C7562-150	Reagent	10 x 15 ml

RELATED PRODUCTS

C7562-12	CK-MB Controls
----------	----------------