



CRP (hs) Wide Range

METHOD

Immunoturbidimetric

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	0.1-3.20.0 mg/dL
Wavelength:	570 nm
Expected Values:	3.0 mg/dL

ORDERING INFORMATION

CATALOG #		SIZE
C7568-40	R1 Buffer	1 x 20 ml
	R2 Antibody	1 x 20 ml
C7568-100	R1 Buffer	1 x 50 ml
	R2 Antibody	1 x 50 ml

RELATED PRODUCTS

C7568-CAL CRP	Calibrator Set
C7568-CTL CRP	Control Set