



Glucose Liquid

METHOD

Hexokinase/G6PD, Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	600 mg/dL
Wavelength:	340 nm
Expected Values:	70-105 mg/dL

ORDERING INFORMATION

CATALOG #	SIZE
G7517-120	1 x 120 ml
G7517-500	1 x 500 ml
G7517-1L	1 x 1 L
Calibrator/Standard Required	

RELATED PRODUCTS

G7518-STD	Glucose Standard
C7506-50	Chemistry Calibrator