



## Glucose Liquid

### METHOD

Oxidase/Trinder, Endpoint

### TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	500 mg/dL
Wavelength:	500 nm
Expected Values:	70-105 mg/dL

### ORDERING INFORMATION

CATALOG #	SIZE
G7521-120	1 x 120 ml
G7521-500	1 x 500 ml
G7521-1L	1 x 1 L
Calibrator/Standard Required	

### RELATED PRODUCTS

G7518-STD	Glucose Standard
C7506-50	Chemistry Calibrator