



Glycohemoglobin

METHOD

Cation-Exchange Resin

TECHNICAL INFORMATION

Storage Temperature:	RT
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	4.0-20%
Wavelength:	415 nm
Expected Values:	6.0-8.3%

ORDERING INFORMATION

CATALOG #		SIZE
G7540-40	Resin	40 test
	Lysing	1 x 120 ml
	Serum Separators	1 x 30 ml
		40
G7540-80	Resin	80 test
	Lysin	1 x 120 ml
	Serum Separators	1 x 30 ml
		80
Standard Required		

RELATED PRODUCTS

G7540-STD	Glycohemoglobin Standard
G7540-2	Glycohemoglobin Controls