



Y-Glutamyl Transferase Liquid

METHOD

Modified Szasz, Kinetic

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	0-800 U/L
Wavelength:	405 nm
Expected Values:	8-54 U/L

ORDERING INFORMATION

CATALOG #	SIZE
G7571-150	1 x 120 ml (R1) 1 x 30 ml (R2)
G7571-450	3 x 120 ml (R1) 1 x 90 ml (R2)
G7571-625	1 x 500 ml (R1) 1 x 125 ml (R2)