



# Hemoglobin A1c

## METHOD

Immunoturbidmetric

## TECHNICAL INFORMATION

Storage Temperature:	02-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	2.0-16.0% mg/dl
Wavelength:	660 nm
Expected Values:	<6.0% non diabetic

## ORDERING INFORMATION

CATALOG #		SIZE
H7541-40	R1 Latex	1 x 30 ml
	R2a Buffer	1 x 9.5 ml
	R2b Antibody	1 x 0.5 ml
	Lyse Reagent	1 x 125 ml
H7541-120	R1 Latex	1 x 90 ml
	R2a Buffer	1 x 28.5 ml
	R2b Antibody	1 x 1.5 ml

## RELATED PRODUCTS

H7541-CAL	HbA1c Calibrator set	
H7541-CTL	HbA1c Control set	
H7541-LYS	HbA1c Lyse Reagent	500 ml