



autoHDL Cholesterol Liquid

METHOD

Direct/Automated

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	150 mg/dL
Wavelength:	600/700 nm
Expected Values:	30-85 mg/dL

ORDERING INFORMATION

CATALOG #	SIZE
H7545-40	1 x 30 ml (R1) 1 x 10 ml (R2)
H7545-80	1 x 60 ml (R1) 1 x 20 ml (R2)
H7545-320	1 x 240 ml (R1) 1 x 80 ml (R2)
H7545-1000	1 x 750 ml (R1) 1 x 250 ml (R2)
Calibrator Required	

RELATED PRODUCTS

H7545-CAL	autoHDL/LDL Calibrator
L7580-18	Lipid Control Set