



# Iron/TIBC

## METHOD

Ferrozine,Endpoint

## TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	500 ug/dL
Wavelength:	560 nm
Expected Values:	
Iron, Total:	60-150 ug/dL
TIBC:	250-400 ug/dL
Iron Saturation:	20-55%

## ORDERING INFORMATION

CATALOG #	SIZE
I7504-60 Iron Buffer UIBC Buffer Iron Color Iron Standard	60 test kit 2 x 75 ml 1 x 120 ml 1 x 15 ml 2 x 30 ml
I7504-400 Iron Buffer UIBC Buffer Iron Color Iron Standard	400 test kit 1 x 1 L 1 x 1 L 1 x 60 ml 1 x 500 ml