



Iron, Total

METHOD

Ferrozine, Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	500 ug/dL
Wavelength:	560 nm
Expected Values:	60-150 ug/dL

ORDERING INFORMATION

CATALOG #		SIZE
17505-120	Buffer Color	1 x 120 ml 1 x 5 ml
17505-500	Buffer Color	1 x 500 ml 1 x 15 ml
<i>Calibrator/Standard Required</i>		

RELATED PRODUCTS

17505-STD	Iron Standard
C7506-50	Chemistry Calibrator