



Phosphorus, Inorganic

METHOD

Phosphomolybdate, Endpoint

TECHNICAL INFORMATION

| | |
|--------------------------|---------------------|
| Storage Temperature: | 2-8°C |
| Reconstituted Stability: | Ready-to-use Liquid |
| Linearity: | 12 mg/dL |
| Wavelength: | 340 nm |
| Expected Values: | 2.5-4.8 mg/dL |

ORDERING INFORMATION

| CATALOG # | SIZE |
|-------------------------------------|------------|
| P7516-120 | 1 x 120 ml |
| P7516-500 | 1 x 500 ml |
| P7516-1L | 1 x 1 L |
| <i>Calibrator/Standard Required</i> | |

RELATED PRODUCTS

| | |
|-----------|----------------------|
| P7516-STD | Phosphorus Standard |
| C7506-50 | Chemistry Calibrator |