



# Microprotein

## METHOD

Pyrogallol Red,Endpoint

## TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	250 mg/dL
Wavelength:	600 nm
Expected Values:	28-141 mg (24 hours)

## ORDERING INFORMATION

CATALOG #		SIZE
P7582-240	Reagent	2 x 120 ml
	Standard	1 x 15 ml

## RELATED PRODUCTS

P7582-STD	Microprotein Standard
P7582-CTL	Urine Chemistry Controls