



Uric Acid

METHOD

Uricase/Trinder,Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	48 Hours RT
	31 Days 2-8°C
Linearity:	25 mg/dL
Wavelength:	520 nm
Expected Values:	2.5-7.7 mg/dL

ORDERING INFORMATION

CATALOG #	SIZE
U7580-150	10 x 15 ml
U7580-300	20 x 15 ml
U7580-400	8 x 50 ml
U7580-1L	4 x 250 ml

Calibrator/Standard Required

RELATED PRODUCTS

U7580-STD	Uric Acid Standard
C7506-50	Chemistry Calibrator