



Uric Acid Liquid

METHOD

Uricase/Trinder, Endpoint

TECHNICAL INFORMATION

Storage Temperature:	2-8°C
Reconstituted Stability:	Ready-to-use Liquid
Linearity:	20 mg/dL
Wavelength:	520 nm
Expected Values:	2.5-7.7 mg/dL

ORDERING INFORMATION

CATALOG #	SIZE
U7581-120	1 x 120 ml
U7581-500	1 x 500 ml
U7581-1L	1 x 1 L
<i>Calibrator/Standard Required</i>	

RELATED PRODUCTS

U7580-STD	Uric Acid Standard
C7506-50	Chemistry Calibrator