

# G6PD REAGENT

## ASSAY FEATURES:

- FDA Cleared / CE Marked
- Kinetic
- Linear to 21.0 U/g Hb
- Correlation  $r = 0.994$
- Stable for 18 months at (2 - 8° C)

A simple, sensitive, and accurate quantitative assay for the kinetic determination of glucose-6-phosphate dehydrogenase (G6PD) in blood.

For use on both automated and manual chemistry analyzers.

## Instrument Applications Available

In the USA - 800.445.9853

Outside the USA - 011.734.487.8300

Email - [marketing@medtest.com](mailto:marketing@medtest.com)

